


# iontoPatch<sup>STAT</sup>®

Electronic Transdermal Iontophoretic Drug Delivery System

The *STAT* means the same IontoPatch benefits in shorter wear-time and...



**More patient compliance** — for those patients who require or prefer shorter wear-time.

**More patient convenience** — can deliver treatment in 4 hours average patient wear-time.\*

(\*Delivery time is dependent upon individual patient skin resistance.)

## Time-released iontophoresis

IontoPatch<sup>STAT</sup>® is a time-released electronic transdermal drug delivery system. An innovative, self-contained battery produces an electric current to carry drug molecules non-invasively across the skin and to underlying tissue. Drug delivery begins upon application and is shut off automatically when the prescribed dosage has been administered. IontoPatch<sup>STAT</sup>® is single-use and disposable.

## Time and cost savings

- Less time to administer in the clinic allowing for other modalities and procedures. Allows more patients to be treated.
- Patients are not tied to external batteries or wires. They can resume daily activities.

## Safety

- Low current and no pH changes reduce skin irritation or burning. There is no need for buffers with competing ions.

## Effective Treatment

- **More** drug delivered — can potentially deliver higher levels of medication over hours than can be delivered in limited time with patient in clinic given patient skin tolerance.
- **More** drug remains localized — low current time-released delivery reduces blood flow to site and keeps drug localized.<sup>1</sup>
- **More** patients improve — low current produces a significantly longer duration of physiologic response.<sup>2</sup>

## Clinician controlled and efficient

- The clinician prescribes, prepares and administers IontoPatch<sup>STAT</sup>®.
- The clinician has more time to provide the patient with additional therapies during the clinic visit.
- IontoPatch<sup>STAT</sup>® is eligible for billing using the existing iontophoretic code.
- No need for or maintenance of external batteries or lead wires.

<sup>1</sup> Sembrowich W et al., *The Drug Delivery Companies Report, Pharma Ventures Autumn/Winter, 53-55, 2003.*

<sup>2</sup> Anderson CR et al., *Physical Therapy, 83: (2), 161-170, 2003.*

Iontophoresis with the battery built in.



Item No.	Product Description	
026501	Original (6) IontoPatch® 80 (80mA-min) w/ (6) saline vials	\$42/box of 6
026502	(6) IontoPatch® SP (40mA-min) w/ (6) saline vials	\$42/box of 6
026503	<b>NEW!</b> (6) IontoPatch STAT® (80mA-min) w/ (6) saline vials	\$48/box of 6